



Presents

An advanced wound care solution,  
Collagen Pad dressing made from  
a native **100% Bovine Type 1 collagen.**



COL-FID4  
CARE™

The power of collagen in  
wound management



**fidia**  
COL-FID4  
CARE™  
Collagen Wound Dressing

INDICATIONS FOR USE: COL-FID 4 Care Collagen Dressing may be used in the management of partial and full thickness wounds, pressure (stage I-IV) and venous ulcers, ulcers caused by mixed vascular etiologies, venous stasis and diabetic ulcers, cuts, and abrasions.  
PRECAUTIONS: • Do not use on people sensitive to bovine products. • Discontinue use and contact a physician if excessive redness, pain, swelling, blistering, infection, no sign of healing or other unexpected symptoms occur. • Avoid eye contact. Flush thoroughly with water if this occurs. • Single use only. Discard any open or unused product. • Store in a cool, dry storage environment, <23°C (77°F).  
DIRECTIONS FOR USE: • Cleanse wound site with a gentle cleanser; apply collagen wound dressing; may be trimmed to size. • Cover with non-adherent dressing; repeat process as indicated. • Dressings should always be changed within 7 days. Frequency of application will depend on drainage. • Collagen wound dressing is intended for short-term use up to 4 weeks. For long-term use, an assessment by a healthcare practitioner is recommended.

1 Dressing | 2x2 in.  
5 Dressings | 2x2 in.

Manufactured for:  
Fidia Pharma USA Inc.  
Farham Park, NJ 07922  
For product inquiries, please call  
1-866-749-2542.

# COL-FID<sup>4</sup> CARE™

The power of collagen in  
wound management

**A-6021**  
(Reimbursement  
code)



**Native 100%  
Bovine Type 1  
Collagen**

**01**

Maintains a moist environment<sup>1</sup>

**02**

Closely resembles the body's native collagen<sup>1</sup>

**03**

Encourages the deposition and organization of new collagen in the wound bed<sup>1</sup>

**04**

Provides a supportive environment for natural wound healing<sup>1</sup>

# COL-FID4 CARE™

The power of collagen in wound management

## SUPPORTS WOUND HEALING

Sourced and manufactured from native 100% Type 1 Bovine collagen, no fillers

Highly absorbent, High-quality collagen

100%



Covers the wound to form a physical barrier and protects it from external environment



Collagen provides a support structure to support the body's natural healing processes.

Easy to apply, non-adherent sterile pad

The collagen pad should be changed within 7 days

For short-term use up to 4 weeks

Each pad is intended for single use

# COL-FID<sup>4</sup> CARE™

The power of collagen in  
wound management

## INDICATIONS FOR USE<sup>2</sup>

Partial and  
full thickness  
wounds



Pressure ulcers  
(stage I - IV)

Venous ulcers



Ulcers caused by  
mixed vascular  
etiologies



Diabetic ulcers



Abrasions  
and cuts




# COL-FID<sup>4</sup> CARE™

The power of collagen in wound management

## INSTRUCTIONS FOR APPLICATION<sup>2</sup>

1

Cleanse wound site with gentle cleanser; apply collagen pad dressing, may be trimmed to size.



2

Cover with non-adherent dressing; repeat process as indicated. Dressings should always be changed within 7 days. Frequency of application will depend on drainage.

3

Collagen dressing is intended for short-term use up to 4 weeks. For long-term use, an assessment by a healthcare practitioner is recommended.



**Once applied to the wound, it transforms into a gel consistency that's adaptable and easily conforms to the wound bed.**

- Absorbs a high amount of exudate and maintains a moist wound environment to aid in the management of wound healing.



# COL-FID4 CARE™

The power of collagen in  
wound management

**A-6021**  
(Reimbursement  
code)



**Native 100%  
Bovine Type 1  
Collagen**

COL-FID4 CARE™ Collagen Dressing may be used in the management of partial and full thickness wounds, pressure (stage I-IV) and venous ulcers, ulcers caused by mixed vascular etiologies, venous stasis and diabetic ulcers, cuts, and abrasions.

**Rx only.**

- Precautions:**
- Do not use on people sensitive to bovine products.
  - This product has not been studied in pediatric populations and is intended only for use in adult patients.
  - This product is not intended for use in procedures where product comes in contact with cerebrospinal fluid.
  - Single use product

**References:**

1. Data on File
2. COLFID4 CARE™ Package Insert



©2025 FIDIA PHARMA USA INC., Florham Park, NJ 07932, a wholly owned subsidiary of Fidia Farmaceutici S.p.A.  
COL-FID4 CARE™ is trademark of Fidia Farmaceutici S.p.A., Abano Terme (PD) Italy  
FID1117-11.2025 November 2025  
For more information, visit our website [www.fidiapharma.us](http://www.fidiapharma.us)